



THERMO MIXER



Thermo Mixer microprocessor based with RS232 interface

Specifications:

Temp Range	Ambient to 99°C
Temperature accuracy	0.5°C
Displacement	6mm
RPM Range	300 – 1500
Thermblock Size	100 × 70mm
Temperature increment	@ 5°C per min
Timer	1 minute to 99 hours 59 minutes, or continuous

Applications: DNA ligation, enzyme reactions or DNA hybridizations etc.

Bio-Age Equipment & Services
S.C.O.-1, Indl. Area Phase-9
Mohali (Pb.)-160062 India
Ph.No.0172-6534210, Fax:0172-5095350
Mail:sales@bioageindia.com, info@bioageindia.com
www.bioageindia.com